
**18<sup>th</sup> INTERNATIONAL BIOLOGY OLYMPIAD**  
**第十八屆國際生物奧林匹亞競賽**  
**JULY 15 - 22, 2007**

*International Biology Olympiad*

*Saskatoon Canada 2007*

**PRACTICAL EXAMINATION 3 實驗題 3**

**Cell Biology/Biochemistry 細胞學/生物化學**

- TASK A. Thiocyanate analysis in cauliflower 29 marks**  
分析花椰菜中硫氰鹽含量 29 分
- TASK B. Determination of the amount of cauliflower needed to be consumed to cause toxicity 5 marks 5 分**  
計算花椰菜中硫氰鹽的毒性
- TASK C. Regulation of gene expression 10 marks**  
基因表現調控 10 分

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**WRITE ALL ANSWERS IN THIS EXAM BOOKLET**  
將答案書寫於本試卷上

**WRITE YOUR 4-DIGIT STUDENT CODE IN THE BOX BELOW AND ON THE TOP OF EACH PAGE OF THIS BOOKLET**  
將四位數的學生代碼填於下面的空格中, 並於每頁的上方重複填寫

<b>STUDENT CODE</b>	
---------------------	--

**Introduction**

The cabbage family contains a class of compounds known as glucosinolates. Some glucosinolates such as glucoraphanin have desired medicinal properties helping to prevent cancers while others such as glucosinalbin have toxic metabolites.

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One of the products of the toxic glucosinolates is the thiocyanate ion ( $\text{SCN}^-$ ).  $\text{SCN}^-$  interferes with iodine metabolism resulting in thyroid hormone deficiency. Eating plants of the crucifer family such as cauliflower will result in the production of a limited amount of thiocyanate ion from glucosinolates such as glucosinalbin.

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The glucosinolate glucoraphanin is metabolized to sulforaphane. Sulforaphane is an inducer of phase 2 proteins. One consequence of phase 2 protein induction is an increased ability of cells to scavenge free radicals and other oxidants. A consequence of decreased oxidant levels is a lower probability of activation of pathways that lead to inflammation. One such pathway is through activation of a protein complex such as NFkappaB.

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**TASK A. To determine the amount of thiocyanate ion released from cauliflower using a spectrophotometric assay. This has a value of 27 marks.**

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**OBJECTIVE:** To use a spectrophotometer to determine how much thiocyanate ion is released from cauliflower. This assay is based upon the principle that in an acid environment thiocyanate reacts with  $\text{Fe}^{3+}$  to form a stable  $\text{Fe}^{2+}$ -SCN red-coloured complex with a maximum absorption at 447 nm.

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**CAUTION: The ferric nitrate reagent solution you will be using is dissolved in 1.0 M nitric acid. Wear gloves and use goggles to protect your eyes before starting the experiment.**

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- Thiocyanate standards in tubes at the following concentrations: 0 micromoles/mL (this is your blank), 0.1 micromoles/mL, 0.5 micromoles/mL, 1.0 micromoles/mL, 2.0 micromoles/mL and 4.0 micromoles/mL.

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- One tube of filtered cauliflower homogenate. 1.0 g of cauliflower was homogenized and the homogenate was diluted to a total volume of 4.0 mL water. This is your unknown and you will be required to determine how many micromoles of thiocyanate are present in one millilitre of this homogenate.  
試管中裝有花椰菜均質液。將 1.0 g 的花椰菜均質化後，最後將均質液的總體積用水調整至 4.0 mL，此即為未知物（為本次實驗所要測定的溶液）。最後要測定的是 1 mL 溶液中含有多少 mM 的硫氰鹽。
- A marker pen to label the cuvettes.  
以馬克筆在分光光度管註記。
- Gloves and protective glasses  
戴上手套與保護鏡。
- On your bench is a spectrophotometer set to an absorbance of 447 nm.  
請注意你桌上的分光光度計是否是設定在 447 nm。

**NOTE:** Before beginning this task, be sure that you have all the materials listed above. If you do not, notify a lab assistant by raising your hand.

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
## Procedure

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1. Put on the gloves and the protective glasses.  
戴上手套與保護鏡
2. To each of the cuvettes containing the ferric nitrate reagent add 100 microliters of each of the thiocyanate standards. The standards are: 0, 0.5, 1.0, 2.0 and 4.0 micromoles thiocyanate/mL. A coloured reaction should become visible except for the 0 micromole thiocyanate standard which serves as your blank. Be sure to label the cuvettes on the frosted surface.  
分別將濃度為 0, 0.5, 1.0, 2.0 與 4.0 M/mL 的硫氰鹽標準液各取 100ul，分別加入分光光度管中。此時，除了濃度為 0 的標準液外，其他會有呈色反應在分光光度管中進行。濃度為 0 的標準液為本實驗的空白對照組。請注意是否有在分光光度管磨砂面上做上註記。

3. To each of the remaining 3 cuvettes add 100 microlitres of the cauliflower homogenate.  
在剩下的 3 個分光光度管中分別加入各 100  $\mu\text{L}$  的花椰菜均質液。
4. Carefully carry the cuvettes to the spectrophotometer which has been set to absorb at 447 nm. Open the lid to the light path in the spectrophotometer and insert the 0 micromole thiocyanate/mL standard (i.e., blank) cuvette. The arrow indicates the light path. Ensure that the walls of the cuvettes through which the light passes is transparent. Close the lid and push the “set reference” button on the top right hand of the panel on the spectrophotometer – see the diagram below.  
**Do not touch any of the other buttons!**
4. 小心地移動分光光度管到分光光度計中，分光光度計的吸收值應該設在 447nm。打開蓋子，插入 0  $\mu\text{M}/\text{mL}$  標準液的分光光度管（本管為空白對照組）。分光光度計裡的箭號，註記的是光束的路徑。請注意光束通過的分光光度管應該是透明的。關上蓋子，按下位於分光光度計操作面板右上方“set reference”按鍵，如下圖所示的位置。  
不可觸碰其他的按鈕。

5. Insert each of the standards and record the reading. Then insert each of the cuvettes containing the unknown and record the spectrophotometer reading. Leave the cuvettes at the spectrophotometer and the laboratory assistants will take care of them.  
分別依序插入其他的標準液的分光光度管，並分別依序記錄下讀數。此時插入未知濃度的花椰菜均質液之分光光度管，並分別記錄讀數。當所有的分光光度管都分別記錄完成後，將分光光度管留置於該處，監考人員會收拾善後。

**Spectrophotometer reading (absorbance) for each standard:** This has a value of 10 marks.

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2.0 micromole/mL thiocyanate: \_\_\_\_\_

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4.0 micromole/mL thiocyanate: \_\_\_\_\_

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**Spectrophotometer reading (absorbance) for the unknown:** This has a value of 3 marks.

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1. \_\_\_\_\_ 2. \_\_\_\_\_ 3. \_\_\_\_\_

6. Plot, on the graph paper (on the next page), the absorbance measurements for your thiocyanate standards against the concentration (micromoles/mL) of the standards. **This has a value of 6 marks.**

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7. Take the average absorbance of your cauliflower homogenate and determine the thiocyanate ion concentration using the previously plotted graph. **This has value of 5 marks.**

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ANSWER: \_\_\_\_\_

8. What is the concentration of thiocyanate in cauliflower? Be sure to state the units (3 marks)

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
ANSWER: \_\_\_\_\_

9. What is the standard deviation of the absorbance of the unknown? **This has value of 2 marks.**

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ANSWER: \_\_\_\_\_

**TASK B. To determine the amount of cauliflower needed to be consumed for it to cause toxic effects because of the presence of thiocyanate (5 marks)**

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**Introduction**

The LD<sub>50</sub> is a toxicology term that describes the dose (i.e. millimoles of toxin/kg animal) of a compound that will kill 50% of the animals tested. In the rat, the LD<sub>50</sub> of sodium thiocyanate consumed is reported to be 9 millimoles/kg. Using the data of the experiment you have just performed, calculate how much cauliflower a rat that weighs 500 g would have to eat in a short time to reach the LD<sub>50</sub> of thiocyanate.

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**Procedure**

Circle the letter of the range that best fits your calculated value. Show your calculations on this page. Continue on the back of this page if necessary.

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**TASK C. To interpret the regulation of gene expression. (12 marks)**

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NFkappaB is a transcription factor complex comprised of two proteins (p50 and p65) bound to a third protein known as IkappaB that is normally present in the cytoplasm. Activation of NFkappaB involves the degradation of IkappaB resulting in the NFkappaB p50/p65 heterodimer translocating to the nucleus where it binds to specific promoter elements increasing the transcription of pro-inflammatory genes such as inducible nitric oxide synthase (iNOS). One indicator of activation of NFkappaB is that the ratio of the p65 to IkappaB protein increases.

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One of the consequences of increased iNOS activity is increased production of the nitric oxide free radical (NO). Nitric oxide reacts with the superoxide anion ( $O_2^-$ ) to form peroxynitrous acid. Peroxynitrous acid is a very strong oxidant.

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Increased oxidant levels often results in activation of NFkappaB while lowering oxidant levels often results in decreased activation of NFkappaB and, hence, lowered levels of expression of pro-inflammatory genes.

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**Procedure**

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1. Examine the figures provided in each of the following sections.

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2. Using the data presented, identify which data set is derived from animals fed a diet high in glucoraphanin and provide the basis for your answer.

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**SECTION A. (5 marks)**

Below is a figure that gives data on NFkappaB activation in spontaneously hypertensive stroke-prone (SHRsp) male rats that were fed one of two diets: a control diet or an experimental diet containing glucoraphanin. In the experimental diet, the animals consumed 10 micromoles glucoraphanin/kg body weight.

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
After several months on these diets, the animals were euthanized, nuclei from the kidneys were isolated and prepared for SDS polyacrylamide electrophoresis. Following separation of the proteins on the gel, the proteins were transferred to nitrocellulose membrane and probed with an antibody that recognized the NFkappaB p65 protein.

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A representative Western blot is shown below (on the left) and next to it is a graph that depicts the quantification of blots from 5 different animals per diet group.

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Answer the following questions:

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1. Which group of animals (**A or B**), were fed the glucoraphanin-containing diet? This has a value of 1 mark.

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ANSWER: \_\_\_\_\_

2. Which of the following statements gives the best explanation for your answer? Circle the letter of that statement. This has a value of 4 mark.

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
**SECTION B.** (5 marks)

Below is a figure that gives data on a 45 kD nitrosylated protein (N-protein) in the kidneys of male SHRsp rats that were put on one of three different diets: a diet containing glucoraphanin and two different control diets.

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The top part of the figure is a representative Western blot while the bottom part of the figure is the quantification of Western blots from 5 different animals per diet group.

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Answer the following questions:

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1. Which of the groups, A, B or C, represent the animals fed a diet containing glucoraphanin? This has a value of 1 mark.

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ANSWER: \_\_\_\_\_

2. Circle the letter of the statement below that best explains your answer. Glucoraphanin treatment results in: (4 marks).

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- THE END -

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**HAVE YOU WRITTEN YOUR STUDENT CODE ON THE FIRST PAGE OF THIS EXAM BOOKLET AND ON THE TOP OF THE OTHER PAGES?**

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